V.Mueller Products and Services

- ENG Cosgrove Flex Clamp
- F Clamp flexible Cosgrove
- D Cosgrove flexible Klemme
- I Pinza flessibile Cosgrove
- E Pinza flexible Cosgrove
- P Pinça flexível Cosgrove
- NL Cosgrove flexibele klem
- S Cosgrove Flex-klämman
- DK Cosgrove Flexibel-klemme

©Copyright 2003, Cardinal Health.
All rights reserved.

V. Mueller, Instra-Lube, and Fogarty-Hydragrip are registered trademarks of Cardinal Health, Inc. or one of its subsidiaries.

Cardinal Health

McGaw Park, IL 60085 USA

Authorized Representative:
Cardinal Health France 205 S.A.S.
F-44110 Châteaubriant

36-2751D

©Copyright 2003, Cardinal Health.

CardinalHealth

Description

The Cosgrove Flex Clamp features a flexible shaft that allows the user a clear field of vision to the surgical site. Additionally, this unique tip design provides an atraumatic occlusion of the vessel. A variety of Fogarty Clamp Inserts are available for use with this clamp. For assembly instructions, refer to manufacturer's directions for use.

Intended Use

Used to occlude a blood vessel temporarily. Used in pulmonary and gastrointestinal procedures and can be used to clamp over indwelling catheters. Also used in minimally invasive and standard open cardiovascular procedures for temporary occlusion of a blood vessel.

For Vascular Work

- Suitable for both veins and arteries.
- Clamps, including the inserts, permit total occlusion of delicate blood vessels with minimal trauma.
- Do not damage surface of vessels and do not penetrate vessel walls.
- Provides occlusion of atherosclerotic vessels without excessive clamping forces.
- Minimizes intimal damage and fragmentation of atherosclerotic material.
- Can be used to clamp over indwelling catheters.

For Pulmonary & Gastrointestinal Work

- Replaces bulky rubber-shod clamps.
- Provides adequate traction and occlusion without crushing the bowel.

The Cosgrove Flex Clamp is designed and manufactured specifically for the purpose of providing pressure and traction as described above and should not be used on hard structures such as bone or cartilage. The jaw of the device should be maneuvered by hand to the vessel site prior to application of jaw closure. Any use of the instrument for any purpose other than that for which it is intended may result in a damaged or broken instrument.

How Supplied

The Cosgrove Flex Clamp is packaged nonsterile.

Fogarty clamp inserts are available through Edward Life Sciences, 1-800-424-3278, or contact your local Cardinal Health distributor.

Inspection

The Cosgrove Flex Clamp should be inspected carefully prior to each use and during the cleaning procedure for any signs of unusual wear, cracking or corrosion. Inspect internal cable for kinking, fraying or any damage to the cable. Microscopic cracks may be detected by the appearance of rust on the instrument. Should any irregularity be noted, do not use.

V. Mueller recommends that this instrument be inspected annually by a qualified V. Mueller repair specialist. The instrument has been etched with the date of manufacture to assist in the tracking of this annual inspection.

Cleaning and Sterilization**Cleaning**

Immediately after each procedure and before sterilization, the Cosgrove Flex Clamp should be thoroughly cleaned with soap and water and brushed to remove any remaining visible debris.

On the Quick Bend styles, find the flushing port located on the bottom of the lower jaw toward the cable end of the jaw. Fill a syringe with cleaning solution and insert the tip of the syringe into the flushing port. Flush thoroughly, working the solution into the jaws by opening and closing the jaws. Repeat flushing.

During manual, ultrasonic or other automated cleaning process with an enzymatic detergent, the jaws of the instrument should be closed partially to permit easy access of cleaning solution to the cable. (Do not clamp ratchet.)

After thoroughly rinsing, carefully inspect and repeat cleaning sequence if any visible debris remains.

 **Caution:** Lubricate the Cosgrove Flex Clamp after each cleaning. Use an approved instrument lubricant such as Instra-Lube concentrate.

Sterilization

The Cosgrove Flex Clamp may be steam sterilized in a peel pouch by standard hospital steam sterilization cycles. Prior to sterilization, the handle must be completely opened and the jaws held manually closed. The following cycles, as defined in AAMI Good Hospital Practice: Steam Sterilization and Sterility Assurance, (ANSI/AAMI ST-46-1993), have been validated for the sterilization of the Cosgrove Flex Clamp:

- Prevacuum Sterilization Cycle (2.5 PSIA pre-exposure evacuation)
Exposure Temperature: 270°F (132°C) Exposure Time: 4 minutes
- Gravity Sterilization Cycles
Exposure Temperature: 250°F (121°C) Exposure Time: 30 minutes
Exposure Temperature: 270°F (132°C) Exposure Time: 15 minutes

The following additional cycle has also been validated for the sterilization of the Cosgrove Flex Clamp:

- Prevacuum Sterilization Cycles (2.5 PSIA pre-exposure evacuation)
Exposure Temperature: 250°F (121°C) Exposure Time: 15 minutes

Product Guarantee

Every surgical instrument bearing the V. Mueller brand name is guaranteed to be free of function defects in workmanship and materials when used normally for its intended surgical purpose. Any V. Mueller instrument proving to be defective will be replaced or repaired at no charge.

Any type of misuse or abuse of the instrument will render the warranty void. V. Mueller assumes no liability if the instrument is misused.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Reference

Atlas of Surgical Instrument Care ©1990, Allegiance Healthcare Corporation, V. Mueller

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