V. Mueller® Products and Services

- ENG Fogarty Clamp
- F Clamp de Fogarty
- D Fogarty Klemme
- I Pinza Fogarty
- E Pinza Fogarty
- P Pinça Fogarty
- NL Fogarty klem
- S Fogarty-klämman
- DK Fogarty-klemme

ENG

Description

A unique grasping principle is incorporated into these clamps to evenly distribute applied pressure while providing adequate traction through the use of replaceable jaw inserts.

A variety of inserts are available for use with these clamps. Softjaw Inserts, Safejaw Inserts, Hydrjaw Inserts gently surround blood vessels and offer differing degrees of atraumatic occlusion. Tractionjaw Inserts have a rigid surface to provide traction while occluding blood flow.

Intended Use

Used to occlude a blood vessel temporarily. Used in pulmonary and gastrointestinal procedures and can be used to clamp over indwelling catheters.

For Vascular Work

- Suitable for both veins and arteries.
- Clamps, including the inserts, permit total occlusion of delicate blood vessels with minimal trauma.
- Do not damage surface of vessels and do not penetrate vessel walls.
- Provides occlusion of atherosclerotic vessels without excessive clamping forces.
- Minimizes intimal damage and fragmentation of atherosclerotic material.
- Can be used to clamp over indwelling catheters.

For Pulmonary & Gastrointestinal Work

- Replaces bulky rubber-shod clamps.
- Provides adequate traction and occlusion without crushing the bowel.

Fogarty Clamps are designed and manufactured specifically for the purpose of providing pressure and traction as described above and should not be used on hard structures such as bone or cartilage. The jaw of the device should be maneuvered by hand to the vessel site prior to application of jaw closure. Any use of the instrument for any purpose other than that for which it is intended may result in a damaged or broken instrument.

How Supplied

The Fogarty Clamp is packaged nonsterile.

Fogarty clamp inserts are available through Baxter Healthcare Corporation, 1-800-424-3278, or contact your local Baxter distributor.

Inspection

The Fogarty Clamp should be inspected carefully prior to each use and during the cleaning procedure for any signs of unusual wear, cracking or corrosion. Microscopic cracks may be detected by the appearance of rust on the instrument. Should any irregularity be noted, do not use.

Cleaning and Sterilization

Cleaning

Immediately after each procedure and before sterilization, the Fogarty Clamp should be thoroughly cleaned with soap and water and brushed to remove any remaining visible debris.

After cleaning and rinsing, carefully inspect and repeat cleaning sequence if any visible debris remains.

Sterilization

The Fogarty Clamp may be steam sterilized in a peel pouch by standard hospital steam sterilization cycles. The following cycles, as defined in AAMI Good Hospital Practice: Steam Sterilization and Sterility Assurance, (SSSA-1988), have been validated for the sterilization of the Fogarty Clamp:

- Prevacuum Sterilization Cycle (2.5 PSIA pre-exposure evacuation)
Exposure Temperature: 270°F (132°C) Exposure Time: 4 minutes
- Gravity Sterilization Cycles
Exposure Temperature: 250°F (121°C) Exposure Time: 30 minutes
Exposure Temperature: 270°F (132°C) Exposure Time: 15 minutes

The following additional cycle has also been validated for the sterilization of the Fogarty Clamp:

- Prevacuum Sterilization Cycles (2.5 PSIA pre-exposure evacuation)
Exposure Temperature: 250°F (121°C) Exposure Time: 15 minutes

Product Guarantee

Every surgical instrument bearing the V. Mueller brand name is guaranteed to be free of function defects in workmanship and materials when used normally for its intended surgical purpose. Any V. Mueller instrument proving to be defective will be replaced or repaired at no charge.

Any type of misuse or abuse of the instrument will render the warranty void. V. Mueller assumes no liability if the instrument is misused.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Reference

Atlas of Surgical Instrument Care ©1990, Allegiance Healthcare Corporation, V. Mueller

©Copyright 2004, 1998, 1993, Cardinal Health.
All rights reserved.

V. Mueller and Fogarty-Hydragrip are registered trademarks of Cardinal Health or one of its subsidiaries.

Cardinal Health

McGaw Park, IL 60085 USA

Authorized Representative:
Cardinal Health France 205 S.A.S.
F-44110 Châteaubriant

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